

DOLOMITE | SYMPHONY



Extra light. Extra adaptable

There is sometimes an obvious need for an extra light walker such as when the user is forced to lift the walker on a daily basis, is weak or has problems breathing.

DOLOMITE SYMPHONY has a rustproof, elegant frame in aluminium, and a total weight of only 6.1-6.8 kg/13.6-15.1 lbs.

The result is an exceptionally adaptable walker that fits just as well in the lounge as it does out in the sea winds.





Aluminium frame

Provides a low product weight and makes it easy to move around with the walker. Its rustproof properties also ensure that the walker remains good looking come rain or shine.

Walking area

Large walking area inside the rollator frame is a signum for **DOLOMITE**. The flip-up seat gives premium walking space.

Makes it easier to reach and access things. The rear forks gives a smooth inside, being friendly to toes and shoes. Structured seat with anti-slip function, rounded rear edge and generous sitting area, means that the user can sit securely and comfortably for a long time.

User friendly brakes

Long brake handle and brake stopper provides good leverage, simple to brake even with weak hand strength. Short distance to park-position reduces risk of losing balance. Short distance to park-

position - reduces the risk of the user losing parking V-shaped brake stopper. Reinforced cable-housing and stainless steel brake stopper.

Features and Options



Tray

Made of impact and resistant plastic. Four stable feet. Can be moved to a table.
Art no. 70311



Basket

Different baskets are available for different height model.
Art no. 12530-15
SYMPHONY 600, 520
Art no. 12538-15
SYMPHONY 450
Max load capacity 7 kg



Cane holder

Makes the cane easily accessible without getting in the way.
Art no. 12606



Backrest

A comfortable and safe. Support to lean against. Mounted on the handlebars and therefore follows the height adjustment.
Art no. 12600

Features and Options



Kerb-climber

"Pedal" that simplifies passing kerbs or similar. Mounted on the right side.
Art no. 12560-21



Hanging bag

Provides extra room when the basket is quickly filled. Can be combined with oxygen tank holder. Art no. 12523



One-hand brake

Makes it possible to activate both brakes with one hand.
Art no. 12507-R (right).
Art no. 12507-L (left)




Brake booster

Provides extra braking force for weak hands.
Art no. 12651



Slow-down brake

A safety feature preventing the walker from rolling away. Perfect for gait training and one handed solutions. Art no 13041-21 (pair).



HP-handle set

For adaptation for one-handed/stroke patients. State size of rollator when ordering. "DOLOMITE HP".
Art no. 12611-R; HP right, one-hand brake left.
Art no. 12611-L; HP left, one-hand brake right



Oxygen cylinder holder

The oxygen cylinder is centrally positioned on the rollator and thereby provides good maneuverability.
Art no. 12682-25 Fits all cylinders.



Adapt-handle set

Can be adapted to completely optional hand position.
Art no. 12613 (pair). State size of rollator when ordering.




Soft grip

A softer alternative to the anatomical handle.
Art no. 12593 (pair)



Wheel guard

Protects clothing from being dirty from wheel splash.
Art no. 12544 (pair)



Crossbar

Make it possible to handle the rollator with one hand.
Art no. 1531573 Right
Art no. 1536813 Left



Screw instead of knob

Replace height adjustment knob with the screw when adjustment is not desired.
Art no. 60140

Three heights



SYMPHONY 450
Art no. 12162



SYMPHONY 520
Art no. 12161



SYMPHONY 600
Art no. 12160

Technical data



Symphony 450 600 mm



Symphony 520 620 mm



Symphony 600 460 mm



Symphony 450 640-790 mm



Symphony 520 480 mm

Symphony 520 600 mm 645 mm 460 mm 705-850 mm 540 mm

Symphony 600 600 mm 710 mm 460 mm 785-935 mm 625 mm



Symphony 450 420 mm



Symphony 520 200 mm



Symphony 450 7,2 kg



Symphony 450 100 kg

Symphony 520 420 mm 200 mm 7,6 kg 100 kg

Symphony 600 420 mm 200 mm 7,8 kg 100 kg

Frame colours



RAL 9006

Invacare Dolomite AB
Växjövägen 303 75
Box 200, 343 75 Diö
Sweden
Tel. +46 (0)476 535 00
Fax +46 (0)476 535 99
www.dolomite.biz

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Symphony - EU - 08/2011



Yes, you can.®